

It's in Your Eyes Awareness and Control HCI Using EOG Signal

System Usage Guide

✚ Welcome to our system! This guide will help you understand how to use the various features and functionalities available. Let's get started.

- **Splash Screen**

- ✓ Upon launching the system, you will be greeted with a captivating splash screen that serves as an introduction to the main home screen. Enjoy the visual experience as the system prepares to assist you.

- **Main Home Screen**

✚ The main home screen is your gateway to different functionalities. Here's an overview of the available options:

1. Virtual Keyboard

- Tap the virtual keyboard button on the main home screen to open the Virtual Keyboard screen. The keyboard is designed to make writing easy and efficient. It features an assisting auto-dictionary, which helps you write faster. Whenever you click on a character, a list of suggested words containing that character will be displayed. These suggestions are generated intelligently, considering the character's position within the word, whether it's at the beginning, middle, or end.
- The Virtual Keyboard offers both an English and Arabic layout. You can seamlessly switch between the two languages by using the dedicated language switch button. Additionally, both keyboards have a

"Speak" button that allows you to hear what you write, aiding in accessibility.

2. Latest News

- The latest news button on the main home screen takes you to the Latest News screen. Here, you will find a curated selection of the most recent news articles. Stay up-to-date with current events by browsing through the headlines and reading the articles that interest you.

3. Chat Inquiries

- Tap the chat inquiries button on the main home screen to access the Chat Inquiries screen. Its functionality is designed to assist doctors in diagnosing patients. The system will present you with a series of questions that require a "Yes" or "No" answer. By responding to these inquiries, doctors can gather valuable information and provide accurate diagnoses.

4. Daily Activities

- The daily activities button on the main home screen opens the Daily Activities screen. This screen offers a collection of the most common activities that people engage in daily. By selecting and interacting with these activities, you can conveniently track and record your daily routines. Additionally, you have the option to use voice commands to interact with the system and log your activities effortlessly.

• Eye-Tracking Cursor

- ✚ Our system incorporates an innovative Eye-Tracking Cursor that responds to Electrooculography (EOG) signals. By monitoring your

eye movements, the system enables you to navigate through the interface and interact with different elements. Look at the desired item or button on the screen, and the Eye-Tracking Cursor will automatically adjust to your gaze, allowing you to make selections intuitively.

 Enjoy using our system, and feel free to explore the various features and functionalities. If you have any questions or need assistance, please contact our support team. We hope this system enhances your user experience and provides valuable assistance in your daily activities.
